

Multisizer 4e data: C:\cell_counter_results\Felix\JF_PBR_day42_T4_02.#m4
Preference file: C:\Multisizer4e\SOP\Default.prf
File ID: JF_PBR_day42_T4
Comment: 50uL sample
Run number: 828
Electrolyte: BCI ISOTON II
Dispersant: None
Aperture: 30 μ m Kd: 44.324
Aperture current: 600 μ A Preamp gain: 4
Size bins: 400 from 0.6 μ m to 18 μ m, log diameter
Total count: 1491 (Coincidence corrected)
Count > 0.6 μ m: 1535 Coincidence corrected: 1539
Coincidence correction: 0.3%
Control mode: Volumetric, 50 μ L
Elapsed time: 13 seconds
Acquired: 19:40 8 Apr 2019
Electrolyte volume: 10 mL
Analytic volume: 50 μ L
Sample: 0.05 mL

Number Statistics (Arithmetic)

JF_PBR_day42_T4_02.#m4

Calculations from 0.600 μ m to 18.00 μ m

Number: 1491
Mean: 0.893 μ m 95% Conf. Limits: 0.866-0.920 μ m
Median: 0.712 μ m S.D.: 0.53 μ m
Mode: 0.608 μ m

d₁₀: 0.614 μ m d₅₀: 0.712 μ m d₉₀: 1.452 μ m

Number Statistics (Arithmetic)

JF_PBR_day42_T4_02.#m4

Calculations from 0.600 μ m to 18.00 μ m

Number: 5.994e6 per mL
Mean: 0.893 μ m 95% Conf. Limits: 0.866-0.920 μ m
Median: 0.712 μ m S.D.: 0.53 μ m
Mode: 0.608 μ m

d₁₀: 0.614 μ m d₅₀: 0.712 μ m d₉₀: 1.452 μ m

Volume Statistics (Arithmetic)

JF_PBR_day42_T4_02.#m4

Calculations from 0.600 μm to 18.00 μm

Volume: 1647 μm^3
 Mean: 3.259 μm 95% Conf. Limits: 3.152-3.366 μm
 Median: 2.387 μm S.D.: 2.11 μm
 Mode: 6.857 μm

d₁₀: 0.799 μm

d₅₀: 2.387 μm

d₉₀: 6.829 μm

Volume Statistics (Arithmetic)

JF_PBR_day42_T4_02.#m4

Calculations from 0.600 μm to 18.00 μm

Volume: 6.622e6 μm^3 per mL
 Mean: 3.259 μm 95% Conf. Limits: 3.152-3.366 μm
 Median: 2.387 μm S.D.: 2.11 μm
 Mode: 6.857 μm

d₁₀: 0.799 μm

d₅₀: 2.387 μm

d₉₀: 6.829 μm

Differential Volume (Smoothing=3)



