

Multisizer 4e data: C:\cell_counter_results\Felix\JF_PBR_day42_T1_02.#m4
Preference file: C:\Multisizer4e\SOP\Default.prf
File ID: JF_PBR_day42_T1
Comment: 50uL sample
Run number: 822
Electrolyte: BCI ISOTON II
Dispersant: None
Aperture: 30 μ m Kd: 44.324
Aperture current: 600 μ A Preamp gain: 4
Size bins: 400 from 0.6 μ m to 18 μ m, log diameter
Total count: 2274 (Coincidence corrected)
Count > 0.6 μ m: 2313 Coincidence corrected: 2321
Coincidence correction: 0.4%
Control mode: Volumetric, 50 μ L
Elapsed time: 13.33 seconds
Acquired: 16:58 8 Apr 2019
Electrolyte volume: 10 mL
Analytic volume: 50 μ L
Sample: 0.05 mL

Number Statistics (Arithmetic)

JF_PBR_day42_T1_02.#m4

Calculations from 0.600 μ m to 18.00 μ m

Number:	2274		
Mean:	0.761 μ m	95% Conf. Limits:	0.745-0.776 μ m
Median:	0.704 μ m	S.D.:	0.38 μ m
Mode:	0.634 μ m		
d ₁₀ :	0.617 μ m	d ₅₀ :	0.704 μ m
		d ₉₀ :	0.899 μ m

