

Multisizer 4e data: C:\cell_counter_results\Felix\JF_PBR_day42_T7_02.#m4
Preference file: C:\Multisizer4e\SOP\Default.prf
File ID: JF_PBR_day42_T7
Comment: 50uL sample
Run number: 834
Electrolyte: BCI ISOTON II
Dispersant: None
Aperture: 30 μ m Kd: 44.324
Aperture current: 600 μ A Preamp gain: 4
Size bins: 400 from 0.6 μ m to 18 μ m, log diameter
Total count: 7136 (Coincidence corrected)
Count > 0.6 μ m: 7229 Coincidence corrected: 7308
Coincidence correction: 1.1%
Control mode: Volumetric, 50 μ L
Elapsed time: 13.24 seconds
Acquired: 19:47 8 Apr 2019
Electrolyte volume: 10 mL
Analytic volume: 50 μ L
Sample: 0.05 mL

Number Statistics (Arithmetic)

JF_PBR_day42_T7_02.#m4

Calculations from 0.600 μ m to 18.00 μ m

Number: 7136
Mean: 0.754 μ m 95% Conf. Limits: 0.748-0.761 μ m
Median: 0.692 μ m S.D.: 0.27 μ m
Mode: 0.603 μ m

d₁₀: 0.615 μ m d₅₀: 0.692 μ m d₉₀: 0.895 μ m

Number Statistics (Arithmetic)

JF_PBR_day42_T7_02.#m4

Calculations from 0.600 μ m to 18.00 μ m

Number: 28.69e6 per mL
Mean: 0.754 μ m 95% Conf. Limits: 0.748-0.761 μ m
Median: 0.692 μ m S.D.: 0.27 μ m
Mode: 0.603 μ m

d₁₀: 0.615 μ m d₅₀: 0.692 μ m d₉₀: 0.895 μ m

Volume Statistics (Arithmetic)

JF_PBR_day42_T7_02.#m4

Calculations from 0.600 μm to 18.00 μm

Volume: 2772 μm^3
 Mean: 1.852 μm 95% Conf. Limits: 1.822-1.881 μm
 Median: 1.079 μm S.D.: 1.29 μm
 Mode: 4.757 μm

d₁₀: 0.650 μm

d₅₀: 1.079 μm

d₉₀: 3.968 μm

Volume Statistics (Arithmetic)

JF_PBR_day42_T7_02.#m4

Calculations from 0.600 μm to 18.00 μm

Volume: 11.14e6 μm^3 per mL
 Mean: 1.852 μm 95% Conf. Limits: 1.822-1.881 μm
 Median: 1.079 μm S.D.: 1.29 μm
 Mode: 4.757 μm

d₁₀: 0.650 μm

d₅₀: 1.079 μm

d₉₀: 3.968 μm

Differential Number (Smoothing=3)



